lionelo

User manual
Instrukcja obsługi
Bedienungsanleitung
Инструкция по эксплуатации
Manuale d'uso
Manuel de l'Utilisateur
Manual de usuario
Handleiding
Naudojimo instrukcija

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








8

Important information

Please read this manual before you start operating the product to learn its functions and use it as intended. The failure to follow the safety instructions and use in compliance with the manual may result in a bodily injury.

1. This device is not a toy. Do not let children play with it. Keep out of reach of children and pets.
2. Check the device regularly for any signs of damage. Do not use a damaged device. Do not try to repair the device yourself. Contact an authorized repair service of the manufacturer. If the power cord is damaged, it should be replaced by the manufacturer to avoid any hazard.
3. Do not use medical substances or medications that are oily or contain essential oils, including herbal extracts, for nebulization. It is acceptable to use suspended liquids if the parameters of the substance are within the values specified in the product specifications.
4. The type, dose and method of drug administration should be consulted with a doctor.
5. Administering the drug to children and people with special needs must take place under close supervision and according to the doctor's guidelines.
6. The device can be used solely for nebulization. Do not use it for any other purpose than its designed use.
7. Before you use the device, consult a doctor.
8. If you experience any discomfort, stop using the device immediately and contact your doctor.
9. Use the device with medicinal drugs which solve in water or contain alcohol and saline solutions. Using any other drugs may cause bronchospasm.
10. The device shall not be used to anesthetize the respiratory system.
11. Before use or a prolonged period of no use, clean and disinfect the medicine cup and accessories.
12. Stop using the device if any of its components are damaged or immersed in water accidentally.
13. Do not switch the device on unless the mask fits the face correctly.
14. Always use any spare parts supplied or recommended by the manufacturer.
15. The device must be switched off and disconnected from the mains after use and before assembly, dismounting and cleaning.
16. Do not immerse the device or its power cord in water or other liquids.
17. After use, remove any remaining liquid from the device.
39. Don't use the unit near high frequency products or electronic products.
40. The charging interface is only used for power supply and cannot be connected to other devices.
41. Performance information provided by the this user manual in accordance with EN ISO 27427:2023 may not apply to drugs supplied in high viscosity form. In such cases, information should be sought from the drug supplier.
42. The highest rated operating altitude:1300m
43. The materials used in the components could not be compatible with solutions/suspensions/emulsions that have not been evaluated.
44. The medicines that have been evaluated:Ambronol hydrochloride solution for inhalation and Salbutamol sulfate solution for inhalation
18. When using the device, maintain the liquid level ensuring that the mesh disc is covered completely.
19. Do not use the device in the presence of any flammable gases or mixtures of oxygen and anesthetics.
20. Do not use the device at high temperatures.
21. Do not use the nebulizer near any heating device or open flame.
22. Do not expose the device to strong shock and vibrations.
23. Do not use microwave ovens to dry or disinfect the device due to the fire hazard.
24. Do not touch the mesh disc with your hand or any sharp objects due to the risk of damage.
25. Do not store the device in a damp or dusty environment.
26. Do not use the device if it was dropped, exposed to heavy moisture or immersed in water. Contact an authorized repair service of the manufacturer.
27. Do not expose the device and its accessories to corrosive liquids and gases.
28. Do not wrap the power cord round the device.
29. Do not remove the battery yourself.
30. If you do not use the device longer than one month, it should be charged at least once a month.
31. Always use the power cord specified in the operation manual.
32. Charge the device before first use.
33. Do not repair or maintain the device when it is operated.
34. The device is operated by the patient.
35. The device is a BF applied part.
36. When you have any doubts concerning switching the device on, using or maintaining it, and also if there is any failure of the device, contact an authorized repair service of the manufacturer.
37. Do not use the device if it begins to malfunction, overheats excessively, shows discoloration or bulging, emits unusual sounds or smells, or if any other abnormal phenomena occur. In such cases, contact an authorized service center.
38. Do not use or leave the device in direct sunlight.
45. One device suggest for single patient use and can be reusable by a single patient.
46. Applied parts not intended to supply heat to a patient.
47. Incorrect replacement of lithium batteries by inadequately trained person would result in an unacceptable risk and hazardous situation.
48. The small parts of this device may be swallowed by the infants and children.
49. The Type-C cable should be kept away from children or out of the reach of children to avoid strangulation.
50. The accessible materials used in the nebulizer are safe for common normal people. For very few operators with extreme skin sensitivity, if any skin discomfort occurs during the use of the nebulizer, please immediately stop using and seek medical advice.

Counter indications

1. Do not use the device with any medications containing pentamidine.
2. The device must not be used by patients with pulmonary edema.
3. The device must not be used by patients who had pulmonary infarction and suffering from acute asthma.
4. If the patient suffers from diabetes or another chronic disease, contact the doctor before use.
5. Do not use if the patient is allergic to aerosol medication.

The product conforms to the requirements of European Union directions.

According to Directive 2012/19/EU, this product is subject to separate waste collection. Do not dispose the product with municipal waste as it may pose hazard to the environment and human health. The product past its service life should be taken to the site recycling electrical and electronic devices.

Intended use


A portable mesh nebulizer is designed for inhalation therapy in health care facilities or for home health care. The device can be used both for adults and for children and people with limited capacities provided it is used under adult supervision.

Specification

Power supply	Battery 3.7 V DC 360 mAh	Vibration frequency	130 kHz, Deviation $\pm 10\%$	Medicine cup volume	6 ml
Power consumption	< 4.0 W	Particle size	MMAD < 5 μm	Safety class	BF applied part
Nebulization rate	0.15 ml/min–0.90 ml/min	Product dimensions / weight	6.8 x 6.6 x 4.05 cm / 57.5 g	Service life	5 years for the device, 2 years for accessories

The particle size and the aerosol output, aerosol output rate measured according to the provisions of the EN ISO 27427:2023, Annex C (using breathing simulator) and D (using the multistage cascade impactor) separately.

The nebulization rate is measured with saline 0.9% solution at a temperature of $24 \pm 2^\circ\text{C}$ and a humidity of 45–75%. It can be vary with medication and ambient conditions. The test solution is albuterol 0.1 % (M/V) concentration in 0.9 % sodium chloride solution. Test conditions of a temperature is $24 \pm 2^\circ\text{C}$ and a humidity is 45–75%. And the plot of cumulative size distribution of the results is as follows.

Note:

1) The particle size distribution and aerosol output may vary by combination of product, medication and ambient conditions such as temperature, humidity and atmospheric pressure.

2) Using a solution, suspension, or emulsion different from that recommended by the manufacturer, in particular, a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer. See drug supplier's data sheet for further details.

3) The effect that the following disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for paediatric or infant populations.

4) Essential performance: particle size and the aerosol output, aerosol output rate.

5) The proportion of mist particles with a diameter of 1 μm to 5 μm generated by the nebulizer is >50%.

Kit

The kit includes the main unit, medicine cup, USB power cord, adult mask, child mask and mouthpiece.

Part list (Fig. 1)

A. Container components

1. Cover
2. Medicine cup
3. Spray nozzle

4. Main unit

5. Power switch
6. USB port
7. Electrode contacts

8. Child mask

9. Adult mask
10. Headband
11. Cup unlocking button

12. Mouthpiece

13. Mesh disc

Assembly

- Attach the container to the main device unit. The container is attached correctly when you hear a characteristic click (Fig. 2). Next, install the appropriate mask or mouthpiece (Fig. 3).

PLEASE NOTE! Before use, clean, disinfect and dry all the device parts as instructed in the "Cleaning and disinfection" chapter.

Power supply

This device has a Type-C cable for charging and powering. It does not come with a power adapter. Please use an AC adapter (output: DC 5.0V 1.0A) approved by IEC 60601-1 for charging and powering. It supports charging while working. It shall be disconnected from the power source after use.

Battery charging

- The fully-charged battery may supply power to the device for up to 60 minutes continuously.
- Low battery level is indicated by a LED which is going to blink blue 5 times and then switch off.
- Use the charger to charge the battery. Ca. 2 hours are required to charge it completely.
- The LED blinks when charging and is solid when it is fully charged.

Operation manual (Fig. 5)

1. LED

LED is solid blue	Operating	LED blinks blue 5 times	Low battery / Switching off
LED blinks blue 3 times	No liquid / The device is not vertical / Switching off	LED blinks green	Charging
LED blinks blue 10 times	The predefined operating time of 10 minutes has elapsed – the device switches off	LED is solid green	Battery fully charged

2. Preparation: Before use, clean, disinfect and dry all the device components as instructed in the "Cleaning and disinfection" chapter.

3. Pouring liquid: open the cover, inject the solution and close the cover. Please note! Before you switch the device on, pour the liquid into the cup.

Please note: protection from leakage! When you pour the medicine into the medicine cup, ensure that you fill it up to the maximum value (6 ml). The recommended liquid amount is 2 to 6 ml. Nebulization takes place solely when the substance to be produced has contact with the mesh disc. Otherwise, nebulization stops automatically. Keep the device upright.

4. Nebulization

- Before you start nebulizing, shake the device horizontally to stir the solution. Based on your individual needs, use the following two nebulization methods (see Fig. 6).
- Put on the mask or place the mouthpiece between your lips and press the power button.
- A – Nebulizing with a mouthpiece B – Nebulizing with a mask
- Take a slow, deep breath and wait until the medication release starts.
- The nebulizer switches off automatically in 10 minutes. To extend the use, press the power button. Ensure there is sufficient liquid in the medicine cup.
- After nebulization is completed, press the power button to switch the device off. Press the container unlocking button, remove the container and empty it (Fig. 7).
- To install the container, press the "press" button and slide the container in the device vertically. Ensure that the container is located correctly (Fig. 8).
- Please note! The liquid may start thickening round the nozzle and the mesh disc which may affect the nebulization result. Stop nebulizing, remove the mouthpiece and other accessories and then wipe the residues away with a sterile gauze. Do not touch the mesh disc with your hand or sharp objects due to the risk of damage.

Manual Cleaning

Following conditions of manual cleaning can be performed for the product: 0.5%-1% neodisher MediClean forte, 10-30 minutes immersion time at 30°C.

If other cleaners are used, select immersion times, concentrations and temperatures of the solutions as instructed by the manufacturer of the cleaner.

- a. Wear disposable protective gloves.
 - b. Put the device in a suitable cleaning basin.
 - c. Create a detergent solution as instructed by the manufacturer of the detergent.
 - d. Completely immerse the device in the detergent solution. Make sure that there are no air bubbles on the products while immersed.
 - e. Dip a clean soft lint-free cloth in the detergent solution, and gently wipe all surfaces and lumen of the device until the device is visually clean. Pay special attention to the product gaps. When the cloth becomes obviously contaminated, discard it, and distribute a new cloth as needed. Then flush the lumen of the device with the detergent solution by using a syringe. Repeat the step for 5 times. The cleaning process takes about 10 minutes.
 - f. Remove the device from the cleaning basin. Rinse the device under running water for 30 seconds at room temperature.
 - g. Put the cleaned and rinsed device into a suitable clean basin with fresh de-ionized water. Thoroughly rinse and flush the lumen of the device with the water for 5 minutes.
- Dry the device through lint free towel.

Manual Disinfection

Following parameters are used:

- Disinfectant: MetriCide™ OPA Plus;
 - Concentration: Ready-to-use solution, undiluted;
 - Temperature: ambient temperature;
 - Soaking time: ≥ 5 min.
 - a. Immerse the accessories in the disinfectant solution.
 - b. Wipe the surface of the device with a sterile soft lint-free cloth to remove any air bubbles on the surface.
 - c. After completion of holding time, remove the device from the disinfection solution.
 - d. Rinsing the devices after manual disinfection with sterile water.
 - e. Wipe and dry the device surface with a sterile lint-free cloth until the surface is free of water stains.
- The disinfection conditions described above are only recommended parameters, please refer to instructions for use of the disinfectant agent for actual use parameters and precautions.

Drying

- Shake the medicine cup to remove any remaining water from the cup and the mesh disc.
- Leave the accessories to dry for at least 4 hours.
- Ensure that all the accessories are perfectly dry, including but not limited to the mesh disc and medicine cup.

Functional Testing, Maintenance:

- Visual inspection for cleanliness of the products.
- All devices should be checked again for dryness.

- After cleaning and disinfection, a thorough inspection and maintenance ensures that the products are fit for use.
 - Check that the product has no dents, cracks, deformations, scratches, etc.;
 - Check all markings on the product for clear visibility.
- Discard and replace any components as necessary.
- Do not use the device with following defects: material deformation, cracks on the product, brittle or other change in the material, etc.
- Test and maintain the device acc. to manufacturer's instruction.

Functional test:

1. Visually inspect each part of the system for cracks or damage, and replace any defective components.
2. Pour 1-6 mL of normal saline (0.9%) into the nebulizer.

3. Functional test:

- Press the Stand-by button, ensure the indicator is lit, and verify that aerosol is visible.
 - Press the Stand-by button again to turn the system off and confirm the indicator is off.
- Discard any remaining liquid before patient use.

Storage and maintenance

1. Storage conditions

Ambient temperature: -20°C--55°C Relative humidity, non-condensing: 10-93% RH Ambient pressure: 70-106 kPa

After removing the product from a temperature outside the specified range, leave it at room temperature for 2 hours.

Other: non-corrosive gases, good ventilation, avoiding high temperature, humidity and direct sun rays.

2. Storage instructions

The device has a shelf life of 5 years, provided the above storage conditions are ensured. The medicine cup should be replaced after no more than 18 months.

After use, clean and disinfect the device and the accessories. Store the device in the original packaging. Avoid any shock.

3. Operating conditions

Temperature: 5°C-40°C Relative humidity, non-condensing: 15-80% RH Ambient pressure: 86-106 kPa

Power supply: lithium battery 3.7 V DC 360 mAh Power supply parameters (not included): 5V DC 1A (power supply unit complying with IEC 60601-1 and 2MOPP requirements only)

Troubleshooting

Problem	Cause/solution
Do not work when turn on.	• Check if there is enough power. • Check if there is enough medication. • Check if the button is function well. • Clear the clogged medication in the vibrating mesh and restart the power.
Low nebulization rate	• Check if the medication cup been filled with right medication, which should be water-solubility, non- corrosive medication. • Check if the medication cup been filled with right volume. • Tilt Main Unit, so that the medication can contact with the vibrating mesh. • Re-assemble the medication cup correctly and restart the device. • Clean the medication cup and vibrating mesh. If it still cannot be used after cleaning, please check whether the vibrating mesh is broken or not. • Clear the clogged medication in the vibrating mesh and restart the power.
After power on, it power off immediately.	• Re-assemble the medication cup correctly and restart the power. • Tilt Main Unit, so that the medication can contact with the vibrating mesh. • Clear the clogged medication in the vibrating mesh and restart the power.
The mesh is clogged with medication.	• Take off the medicine cup. wash with warm water and air dry.
The medication cup assembly leakages.	• Re-assemble the medication cup correctly. • Contact the manufacturer or the distributor for a new one.
Residual medicine liquid in the medication cup.	• It is a normal phenomenon.Clean the medication cup after using.
Can not charge.	• Make sure the Type-C cable,the adapter are connected well. • Contact the manufacturer or the distributor for a new one.

Specific indications for children

- The mask should cover the mouth and the nose.
- Do not nebulize a sleeping child.
- Please note! Children require adult help and supervision when using the device.

Electromagnetic compatibility (EMC) symbols

1. This device meet the requirement of electromagnetic compatibility in IEC 60601-1-2.
 2. The user needs to install and use according to electromagnetism compatibility information which is attached with it.
 3. Portable and mobile RF communication device and some household appliances, such as mobile, interphone, microwave oven, dry blower, may influence Portable Mesh Nebulizer performance, so Portable Mesh Nebulizer should be kept away from them during using.
 4. Guidance and manufacturer's declaration stated in the appendix.
1. With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. 2. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. 3. Medical devices should also not interfere with other devices. 4. In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the IEC 60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. 5. Medical devices manufactured by FEELLIFE HEALTH INC. conform to this IEC 60601-1-2 standard for both immunity and emissions. Nevertheless, special precautions need to be observed:

Note: The EMISSIONS characteristics of this equipment make it suitable for use in a residential environment (for which CISPR 11 class B is normally required).

Warning:

1. The use of accessories and cables other than those specified by FEELLIFE, with the exception of cables sold by FEELLIFE as replacement parts for internal components, may result in increased emission or decreased immunity of the device.
2. The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.
3. PORTABLE RF communications equipment (including peripherals as such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance this equipment could result.
4. Do not use mobile (cellular) telephones and other devices, which generate strong electrical or electromagnetic fields, near the medical device. This may result in incorrect operation of the unit and create a potentially unsafe situation. Recommendation is to keep a minimum distance of 7 m. Verify correct operation of the device in case the distance is shorter.

Details of cable:

Name	Type	Length(m)	Cable Shielded
Type-C Cable	Type-C	1.0	UnShielded

Guidance and Manufacturer's declaration – electromagnetic emissions		
The device are intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance & Declaration – electromagnetic immunity

The device are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity testing	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line to line ±0.5 kV, ±1kV, ±2kV line to ground	±0.5 kV, ±1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT) for 1 cycle 70% UT (30% dip in UT) for 25/30 cycles <5% UT (>95 % dip in UT) for 5 sec	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT) for 1 cycle 70% UT (30% dip in UT) for 25/30 cycles <5% UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device require continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The magnetic field with the feeder frequency should have the level corresponding to the standard commercial or hospital environment.
Proximity magnetic fields	CW 8A/m for 30kHz Pluse modulation 2.1kHz,65A/m for 134.2kHz Pluse modulation 50kHz,7.5A/m for 13.56MHz	CW 8A/m for 30kHz Pluse modulation 2.1kHz,65A/m for 134.2kHz Pluse modulation 50kHz,7.5A/m for 13.56MHz	Power frequency magnetic fields/Proximity magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance & Declaration – electromagnetic immunity

The device are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such an environment.

Immunity testing	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidelines
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM and amateur radio bands	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM and amateur radio bands	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 & 10 V/m 80 MHz to 2.7 GHz 385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014+A1:2020)	3 & 10 V/m 80 MHz to 2.7 GHz 385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014+A1:2020)	If higher IMMUNITY TEST LEVELS than those specified in Table 9 are used, the minimum separation distance may be lowered. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = [6/d] \times \sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

Software version: WHQA14 V1.1

The photos are illustrative and the actual appearance of the products may differ from the one in the photos.

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EN: The symbol placed on the product indicates that the batteries or accumulators used in it should not be disposed of with household waste. Such waste must be properly sorted in accordance with applicable regulations. Used batteries and accumulators must not be discarded in unsorted waste, as they may pose a risk to the environment and human health. Users should utilize available selective collection points, which enable their return, recycling, and safe disposal. For detailed information on local battery and accumulator recycling regulations, please contact your municipal office, local waste management facility, or landfill.

Note: Lithium batteries may pose a fire or explosion hazard, especially if damaged, overheated, or improperly disposed of. The chemicals contained within them may also be harmful to health and the environment.

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Lotus NL B.V.